



Plates Preparation

Designing Life with AI

July 14, 2026

Description

Method for making nutrient agar plates in petri dishes for bacterial growth.

Material

- LBA medium (LB + agar)
- Antibiotic stock solution (if required; 1:1000 final ratio)
- Petri dishes

Steps

Microwave the medium carefully: slightly loosen the cap of the LBA bottle and heat for 30 s, let it cool briefly, then repeat in 30 s intervals until the LBA is fully melted.

1. Let the agar cool down a bit.
2. If adding antibiotic, put x μ l of antibiotic for x ml of LBA (want a ratio of 1:1000) and label the plates with a unique code on the side.

The following steps have to be performed under highly sterile conditions (under flame or hood) to avoid contamination:

3. Pour the melted agar media in the plates, fill the plates to half.
4. Move the plates in small circles to make the bubbles drift to the side.
5. Let it solidify fully.
6. Store in a plastic bag, labelled with date and potential antibiotic.